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Portable induction loop amplifier

ENGLISH

eLoop / eLoop+

Portable induction loop amplifier

EN

Installation and user manual
eLoop / eLoop+

Table of contents

1. Introduction.....	16
1.1 Purpose	16
1.2 Target audience	16
1.3 Alerts	16
1.4 Icons.....	16
1.4.1 Notes icons.....	16
1.4.2 "Caution", "Warning" and "Danger" icons.....	16
2. Description	17
2.1 eLoop.....	17
2.2 Contents.....	17
2.3 eLoop operation.....	18
2.4 Labeling system.....	19
3. Installation	21
4. Set your eLoop into service	21
4.1 Battery.....	21
4.2 Microphone jack for OP-M80.....	21
4.3 Handset jack plug for OP-E.....	21
4.4 Putting the eLoop out of service.....	21
4.5 Recharging the battery.....	21
4.6 Regulations.....	21
5. Technical specifications	22
6. Informations.....	23
6.1 Maintenance and care.....	23
6.2 Warranty.....	23
6.3 Disposal of used electric and electronic unit	23
6.4 Rechargeable battery.....	23
6.5 CE certification.....	23

1. Introduction

1.1 Purpose

The installation and user manual provides the information needed to install, to configure and to use your eLoop.

1.2 Target audience

The installation and user manual is destined for eLoop installers and users.

1.3 Alerts

This manual mentions four types of warning. These warnings indicate a potentially hazardous situation which, if mishandled, could result in serious, moderate or minor personal injury, and/or property damage

- **Note**

Additional information. Generally, failure to comply with a "Remark" alert will not result in property damage or personal injury.

- **Caution**

Failure to observe a "Caution" alert may result in property damage.

- **Warning**

Failure to comply with a "Warning" type alert may result in serious property damage and personal injury.



- **Danger**

Failure to comply with a "Danger" alert can result in death.

1.4 Icons




1.4.1 «Notes » icons

The "Notes" icons provide additional information about the note. See the following examples :

1.4.2 "Caution", "Warning" and "Danger" icons

The "Caution", "Warning" and "Danger" icons indicate the type of risk present. See the following examples:

2. Description

eLoop is a product designed to ensure audio accessibility and one to one communication for people suffering hearing loss. This product is destined to all establishments receiving public. It has the advantage of being a portable solution that does not require prior installation : it can be directly affixed to the furniture. Easy to use, eLoop ensures an optimal result and a durable installation.

2.1 eLoop

eLoop is a portable induction loop that can be used to equip counters, reception desks, intercom systems, small rooms, etc. The system meets the equal opportunities law restrictions while complying with the requirements of the EN60118-4 standard. The set is designed to offer accessibility for hard of hearing people equipped with a hearing aid with the T position.

The amplifier is designed to be placed directly on a desk or counter. It is fitted with a rechargeable battery, two inputs, settings and LEDs.

The device integrates audio processing to automatically control the recording levels, avoiding loud and sudden noises.

2.2 Contents

According to the version, the kit includes:

- **eLoop version**
 - An eLoop amplifier
 - A rechargeable battery
 - A power supply
 - Set of 2 stickers “Adapted for the hard of hearing people area”.

- **eLoop+ version**
 - An eLoop amplifier
 - A rechargeable battery
 - A power supply
 - Set of 2 stickers “Adapted for the hard of hearing people area”.
 - An OP-E3 handset
 - An OP-M80 gooseneck microphone


Check that all the items listed above are included in the box. If not, please contact us.

2.3 eLoop operation


eLoop is a portable amplifier that integrates a magnetic loop. This system allows direct communication between an interlocutor and an hard of hearing person with an hearing aid equipped with the "T" or "T-coil". Embedded into the hearing aid, the coil consists of a copper coil. Thus, the magnetic field generated by eLoop will pass through the small loops in the hearing aid and, by the induction principle, the electrical signal present in the large loop will find its way into the small loops. The signal is transmitted from the audio amplifier to the hearing aid, which then returns it to the hard of hearing people's ears.

To guarantee an optimal use of the eLoop it is necessary to put it on the furniture between the interlocutor and the user.

The presence of an earring loop is signaled by a blue logo representing a crossed-out ear and a letter T. Generally, an hearing aid has two major positions, the M and the T. The M position allows sound to be perceived through sound pressure like a microphone, while the T position (T for telephone) receives directly the audio signals transmitted by induction via the integrated coil. Some hearing aids combine these two modes of operation via the MT position. It allows the hard of hearing people to perceive both ambient noise and inductive signals.

2.4 Labeling system


(1) OP-E3 Handset

(2) Front cover

(3) OP-E3 Handset Stand

(4) Back cover

(5) Antenna for induction loop


(6) Battery

(7) Space for battery

(8) eLoop base

(9) Battery cover

2.4 Nomenclature

(10) Handset potentiometer

(11) Handset jack plug

(12) Power supply plug

(13) Charging indicator

(14) Low battery indicator

(15) ON/OFF button

(16) Speech diffusion-indicator

(17) Operation indicator

(18) Microphone capsule (internal)

(19) Microphone jack plug

(20) Microphone potentiometer

3. Installation

eLoop is a ready-to-use system, all you have to do is place it at the right distance between you and the person you are speaking to.

It is also possible to screw the fixing slot of the eLoop, onto your furniture thanks to two dedicated place. This solution is recommended for a single spot usage. To take the eLoop out from its slot, just slide it gently forward.

4. Putting eLoop into service

To turn on e-LOOP, please press the button  at the back of the product.

4.1. Battery

To access the eLoop's battery, slide the cover at the back of the product. It is necessary to prime the battery to use eLoop. To do so, it must be plugged in: connect the power supply to the socket indicated by the pictogram .

4.2. Microphone jack for OP-M80

Available as an option, the OP-M80 microphone captures the sound emitted by the agent's voice and transmits it directly to the hard of hearing people's hearing aid or to the eLoop's integrated earpiece if the hard of hearing people do not have an hearing aid.

To use the OP-M80, plug the microphone jack into the socket marked with the pictogram .

To adjust the microphone level, turn the potentiometer indicated by the pictogram .

4.3. Handset jack plug for OP-E3

Available as an option, the OP-E3 earpiece allows hard of hearing people without hearing aids to receive the sound emitted. To benefit from its use, the user must simply pick up the earphone from its support and place it on his ear.

To use this system, simply plug the earphone jack into the socket indicated by the pictogram .

To adjust the volume of the earpiece, please turn the potentiometer indicated by the pictogram .

4.4. Putting the eLoop out of service

To turn off eLoop, please press the button  at the back of the product.

4.5. Recharging the battery

eLoop can work with its integrated rechargeable battery or with a power supply. The eLoop's battery holds an autonomy of 10 hours and recharges completely after 4 hours of charging. To charge eLoop, plug the power supply into the socket indicated by the pictogram .

LEDs at the back of the product indicate when the battery is low  or charging .

4.6 Regulations

For an optimal use and to comply with the regulations in force, it is necessary to indicate the presence of a solution for the hearing impaired by sticking on the stickers provided.

5. Technical specifications

AUDIO INPUT	
Integrated microphone	Capsule microphone
Microphone plug	1x 3.5mm Micro Jack
Phantom power supply	5V, 1mA
AUDIO OUTPUT	
Headphone plug	1x 3.5mm Micro Jack
POWER SUPPLY	
Voltage / Current	5 Vdc / 1,5A
Connector	Micro USB
Battery	Lithium Ion 3,7V 4000mAh
Autonomy	10 hours
Charge time	4 hours
TECHNICAL SPECIFICATIONS	
Current	3,5A rms
Protections	Thermal resistance Overload
Coverage	1,10m
PRODUCT SPECIFICATIONS	
Dimensions	250x185x80mm
Net weight	366g
Packaging	Cardboard
Shipment dimensions	400x250x100
Shipment weight	1Kg
LEDS DISPLAY	
Information	<p>Charging battery (1) Low battery (2) Output current (3) Power supply (4)</p> 

6. Information

6.1. Maintenance and care

The eLoop does not require any maintenance. If the unit becomes dirty, simply wipe it clean with a soft, damp cloth. Never use spirits, thinners or other organic solvents.

Do not set up the unit where it will be exposed to full sunlight for long periods. In addition it must be protected against excessive heat, moisture and severe mechanical shock.

Note: This product is not protected against splash water. Do not place any containers filled with water, such as flower vases, or anything with an open flame, such as a lit candle, on or near the product.

6.2. Warranty

The eLoop is a very reliable product. Should a malfunction occur despite the unit having been set up and operated correctly, please contact your dealer.

The warranty covers the repair of the product and returning it to you free of charge during the warranty period. It is recommended that you send in the product in its original packaging, so keep the packaging for the duration of the warranty period.

The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by people not authorized to do so (destruction of the seal on the unit). Repairs under warranty are only carried out providing you have filled in and returned the enclosed warranty card from the dealer as well as a copy of the sales slip.

Always specify the product number in any event.

6.3. Disposal of used electric and electronic unit

(Applicable in the countries of the European Union and other European countries with a separate collection system).

The symbol on the product or the packaging indicates that this product is not to be handled as ordinary household waste but has to be returned to a collecting point for the recycling of electric and electronic units. You protect health and environment by the correct disposal of this product. Material recycling helps to reduce the consumption of raw material. You will receive further information on the recycling of this product from your local community, your communal disposal company or your local dealer.

6.4. Rechargeable battery

The rechargeable battery that comes with the device can be recycled. Please dispose of this rechargeable battery in receptacles designed for that purpose or return it to a retail outlet. Dispose only of those batteries that have been completely discharged to ensure environmental protection.

6.5. CE certification

This device complies with the following EC directives:

- 2011 / 65 / EC RoHS directive
- 2012 / 19 / EC WEEE directive
- 2014 / 30 / EC EMC directive
- 2014 / 35 / EC Low voltage directive

Compliance with the directives listed above is confirmed by the CE seal on the device. CE compliance declarations are available on the Internet at www.opus-technologies.fr.

Technical specifications subject to change without notice.

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